


# ToF-ACSM Calibrations


Leah Williams, Phil Croteau, Wen Xu

# Setup for Mass Based IE Determination



*Plot Measured Mass vs Input Mass*

# Aerosol Diluter Parts





- HEPA filter, Pall p/n 12144 (with 3/8" FNPT ports)
- Valve Swagelok SS-1RS4
- Capillary tube .032"ID x 1/8" OD x 8" length
- Two 1/4" swage TEEs, two 1/4" x 1/8" swage tube stub reducer (you will probably want one of these to be a bored thru fitting), two 3/8" MNPT x 1/4" swage, one 1/4" union, two 1/4" port connector.
- Some 1/4" st. st. tubing (straight and with 90 degree bends).
- Mixing tube. McMaster 1/4" OD 3530K43 \$159.31  
(<http://www.koflo.com/static-mixers/stratos-tube-mixers.html>) 7" length



# $\text{NH}_4\text{NO}_3$ Solution Concentration for ACSM IE Calibration


Want to minimize multiply charged diameters from DMA

Tests performed with TSI 3076 Atomizer



- Use 0.005M  $\text{NH}_4\text{NO}_3$
- Aim for 1000 p/cc with diluter valve closed

# Setup for calibration



What are we calibrating? Ionization efficiency (IE).

- Want response of instrument to a known amount of  $\text{NO}_3$ , or how many ions (at  $m/z$  30 and 46) do we detect per pg of  $\text{NO}_3$ . (For the AMS, IE is expressed as ions/molecule.)
  - $\text{NH}_4\text{NO}_3$  has CE of 1.
  - $\text{NO}_3$  forms mostly  $m/z$  30 and 46, so easy to count all the ions.
- Use size-selected  $\text{NH}_4\text{NO}_3$  (AN) particles – 300 nm is optimal.
  - Enough mass to get good signal to noise.
  - Larger sizes might not be fully transmitted by lens.
  - Smaller sizes have many more Q2 particles from DMA.
- Count particles with CPC.
  - Use Analog Input 3 on ToF-ACSM for Aquility to record CPC counts.
  - Remember that a CPC can be off by  $\pm 20\%$ !

# IE\_NO3 Calculation

- $$\frac{\text{IE} = \text{ions/s (30 + 46)}}{\text{pg/s (NO}_3\text{)}}$$
- Add the signals at mz 30 and mz 46 = ions/s (30+46).
- $\text{pg/s} = \text{CPCcnts} * 4/3\pi(D_{\text{mob}}/2)^3 * \rho * \text{SF} * 1\text{e-9} * \text{MFA}_{\text{anion}} * \text{Flowrate}$ 
  - $\rho$  = density (1.72 for  $\text{NH}_4\text{NO}_3$ , 1.77 for  $(\text{NH}_4)_2\text{SO}_4$ )
  - SF = shape factor = 0.8 for  $\text{NH}_4\text{NO}_3$ , 0.85 for  $(\text{NH}_4)_2\text{SO}_4$
  - $\text{MFA}_{\text{anion}}$  = mass fraction anion =  $62/(62+18)$  for  $\text{NO}_3$
- Plot ions/s for (30+46) vs pg/s for  $\text{NO}_3$
- Slope = IE\_NO3 (ions/pg)
- Also record AB and flow rate during IE calibration.

# RIE\_NH4 Calculation

- $$\text{RIE\_NH4} = \frac{\text{ions/pg (NH4)}}{\text{ions/pg (30 + 46)}} = \frac{\text{CPCcnts} * 4/3\pi(\text{Dmob}/2)^3 * \rho * \text{SF} * 1\text{e-9} * \text{MFCation} * \text{Flowrate}}{\text{CPCcnts} * 4/3\pi(\text{Dmob}/2)^3 * \rho * \text{SF} * 1\text{e-9} * \text{MFAnion} * \text{Flowrate}}$$
- Note that in the calculation of pg/s, the only difference for NH4 and NO3 is the MFCation (62/(62+18)) or MFAnion (18/(62+18)).
- $$\text{RIE\_NH4} = \frac{\text{NH4 (ion/s)}/18}{\frac{\text{NO3 ion/s}/62}{1.05}} = \frac{\text{NH4 (ions/s)} * 1.05 * 62}{\text{NO3 (ions/s)} * 18}$$
- Plot NH4 (ions/s) \* 1.05 \* 62 vs NO3 (ions/s) \* 18
- slope = RIE\_NH4



# RIE\_SO4 Calculation




- (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> (AS) has a low (and variable) CE. Reference SO<sub>4</sub> to NH<sub>4</sub> in internally mixed particles, then CE does not matter.
- $$\text{RIE\_SO4} = \frac{\text{ions/pg (so4)}}{\text{ions/pg (30+46)}} = \frac{\text{ions/pg(so4)*ions/pg(nh4fromAN)}}{\text{ions/pg (30+46)*ions/pg(nh4fromAS)}} = \text{RIE\_NH4}$$
- Assumption: RIE\_NH<sub>4</sub> is the same in AS and AN
$$= \frac{\text{ions/pg(so4)} * \text{RIE\_NH4}}{\text{ions/pg(nh4fromAS)}}$$
- $\text{SO4 pg/s} = \text{CPCcnts} * \frac{4}{3}\pi(\text{Dmob}/2)^3 * \rho * \text{SF} * 1\text{e-9} * \text{MFA}_{\text{anion}} * \text{Flowrate}$ 
  - $\rho$  = density (1.72 for NH<sub>4</sub>NO<sub>3</sub>, 1.77 for (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>)
  - SF = shape factor = 0.8 for NH<sub>4</sub>NO<sub>3</sub>, 0.85 for (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>
  - MFA<sub>anion</sub> = mass fraction anion = 96/(96+36) for SO<sub>4</sub>
  - MFC<sub>cation</sub> = 36/(96+36) for NH<sub>4</sub> in (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>

## RIE\_SO4 cont.

- $$\text{RIE\_SO4} = \frac{\text{ions/pg (so4)} * \text{RIE\_NH4}}{\text{ions/pg(nh4fromAS)}}$$
- Note that in the calculation of pg/s, the only difference for NH4 and SO4 is the MFCation or MFAAnion.
- $$\text{RIE\_SO4} = \frac{\text{SO4 (ions/s)/96} * \text{RIE\_NH4}}{\text{NH4 (ions/s)/36}} = \frac{\text{SO4 (ions/s)} * 36 * \text{RIE\_NH4}}{\text{NH4(ions/s)} * 96}$$
- Plot  $\text{SO4 (ions/s)} * 36 * \text{RIE\_NH4}$  vs  $\text{NH4 (ions/s)} * 96$
- slope = RIE\_SO4

# Acquility Calibration Interface



Chose ACSM Calibration profile

Follow steps 1 to 3, then hit Start.


After data collected, will ask for CPC counts (if not recorded on AI3), mobility diameter, and species.

Then will ask if want to start a new graph or add to existing graph.

**Timing should be the same as for ACSM\_Aquisition mode.**


- Switch filter/sample valve every 20 sec.
- Save 2 to 5 minute data file.
- Multipoint  $\text{NH}_4\text{NO}_3$  calibration
- Multipoint  $\text{RIE\_SO}_4$  calibration.

# IE\_NO3 and RIE\_NH4 in Acquility



- Do multiple concentrations of NH4NO3, including 0
- IE\_NO3 and RIE\_NH4 are recorded in workspace profile
- Saved in h5 file
- Imported to Tofware

# RIE\_SO4 in Acquility



- Do multiple concentrations of  $(\text{NH}_4)_2\text{SO}_4$ , including 0
- RIE\_SO4 is recorded in workspace profile
- Saved in h5 file
- Imported to Tofware



# Manual IE Calibration

- Highly recommend that users collect calibration data in ACSM\_Acq mode and analyze in Tofware to confirm that get same answer as Acquility auto calibrations.
- Detailed instructions in ToFACSM\_ManualIECalTutorial\_COST.ifn (open notebook in Igor Tofware experiment).
- Acquire data in ACSM\_Aquisition mode, i.e., with filter/ambient switching the same as you will use for data acquisition.
- Take data for several number concentrations of 300 nm  $\text{NH}_4\text{NO}_3$  particles and several concentrations of 300 nm  $(\text{NH}_4)_2\text{SO}_4$  particles, including 0 particles.
- Load into Tofware. Be sure to UNCHECK Correct plots for airbeam (in ACSM Settings).
- Calculate species with  $\mu\text{g}/\text{m}^3$  UNCHECKED. This will give you ions/sec without RIEs applied.
- Equations for IE\_NO3, RIE\_NH4 and RIE\_SO4 are on slides 7, 8 and 10. Also in notebook.

# Notes on cleaning TSI atomizer

- Slow way
  - Use a bottle of pure water, run atomizer until CPC counts get to zero
    - This takes ~10 minutes when you're size and/or mass selecting particles
- Faster way
  - Remove top fitting, flush with water
  - Push water through small tube with a rinse bottle
  - Still need to run pure water after, but not for long.

Tested that RIE\_NH4 same for  $\text{NH}_4\text{NO}_3$  and  $(\text{NH}_4)_2\text{SO}_4$  using mixtures (Xu et al., AST, 2017).



RIE\_NH4 is consistent across range of mixtures, i.e., no matrix effects.